

We are proud to announce that
JMI Infusion Set is now...



LATEX FREE

With

Convenient Y-Injection Port
ISO standard Luer Connector
& Chemical Indicator

... more patient compliant

Dear Valued Customer,

In compliance with the global requirement of quality standard, JMI has updated Infusion Set in order to facilitate ease of extra injection administration, minimization of probable health hazards and more patient compliance.

Introduced Y-Injection Port by replacing Latex tube

- Facilitates ease of extra injection administration
- Minimizes probable health hazards caused by Latex tube
- Ensures more patient compliance
- Establishes safety and durability of the product

Introduction of Luer Slip Connector by replacing conventional one

- Meets Internationally recognized and ISO 594 compliant connection system
- Ensures “leak-proof” connection with all types of Luer fittings e.g. IV Cannula, Butterfly needle, Disposable Needle

... more user friendly







These new features attained trust of the Users

A multi-center clinical trial of JMI Infusion set was conducted in a large number of reputed hospitals of Bangladesh which resulted in an excellent user satisfaction in related parameters.¹

Parameters:

1. Customers appreciated the initiative & planning for upgrading JMI Infusion Set to comply with the global standard requirement
2. Customers appreciated overall development of JMI Infusion Set compared to the conventional ones
3. Customers appreciated overall convenience of IV Fluid administration with upgraded JMI Infusion Set
4. Customers appreciated overall performance of Y-Injection Port of JMI Infusion Set to facilitate ease of extra injection administration
5. Customers appreciated overall safety and convenience of upgraded JMI Infusion Set
6. Customers appreciated “leak-proof” ISO 594 compliant connection system of upgraded JMI Infusion Set

User Satisfaction Level*



* Average

... confirms users satisfaction

Glimpse OF JMI INFUSION SET



Spike with Cap

- Sharp spike for optimum penetration of glass/plastic bottle rubber stopper
- Spike Cap is provided as “safety guard” to prevent unwanted dust penetration
- Available as both Vent and Non-Vent Type
- Delivers 20 Drops / ml (approx.)



Drip Chamber

- Transparent
- Flexible
- Easy to count drops per minute




Tube

- Transparent
- Soft
- Kink resistant



Flow Regulator

- Specially designed roller clamp (flow regulator device) facilitates fine adjustment and control of fluid flow
- Light weight clamp enables single handed operation of the flow regulator device



Y-Injection Port

- Latex Free
- Minimizes probable health hazards (allergic manifestation)
- Facilitates ease of extra injection administration



Luer Slip Connector

- Internationally recognized and ISO 594 compliant
- Ensures “leak-proof” connection



JMI INFUSION SET



With special features:

- LATEX FREE
- Convenient Y-Injection Port
- ISO standard Luer Connector
- Chemical Indicator

Chemical Indicator

... is an unique feature to indicate sterilization status of the products inside Inner Carton and Master Carton of JMI Infusion Set. After EO Sterilization, strip color turns from Pink to Green.

...Confirms safety of the products




Fig.: Chemical Indicator Strip

TOXIKON
ADVANCING YOUR INNOVATION

Moreover,

Bio-compatibility study of JMI Infusion Set was conducted by Toxikon Corporation of USA and the test results were found to be in compliance with the desired standards.²

...Meets the Global Standard



Corporate Office:

7/A, Shantibag, Dhaka-1217, Bangladesh, Tel: +88-02-9346630, 8318375, 9333102, 831733
Fax: +88-02-9337798, 831303, E-mail: info@jmigroup-bd.com, Web: www.jmigroup-bd.com

Ref.: 1.Inhouse Data, 2. Toxikon Corporation Report